

● CLOVER



misano^{MIS}



Complete set of sizes for open, mis and deformity procedures

Self-tapping screw for insertion even without tapping

Great mechanical properties

Double thread for improved sealing

Sophisticated yet compact instrumentation

Cannulated screw for injection of cement or bone substitute

Sterile packaging

Stabilization system thoracolumbar vertebral

Misano is the first stabilization system thoracolumbosacral, totally made of titanium, which allows the use of a single pedicle screw for the treatment of all degenerative, traumatic, and deformity pathologies. The self-tapping screw can be inserted without the aid of the tap, has a specific two-principle thread at the distal level that ensures rapid implantation, which then transitions to four-principles while also ensuring a high pedicle seal. The screws are also always cannulated to allow insertion with the aid of a guide wire and fenestrated to allow cementing.

The Misano thoracolumbosacral stabilization system from Clover Orthopedics consists of single- and polyaxial pedicle screws, clamping nuts, straight and pre-curved titanium and cobalt-chrome bars, a wide range of hooks, offsets, connectors and cross-links, and is used to facilitate a solid arthrodesis of the treated portion of the spine.





Appropriately used, the Misano thoracolumbar-sacral stabilization system from Clover Orthopedics is indicated to promote the development of solid thoracic, lumbar and sacral arthrodesis.

It is recommended in cases of spinal deformity, degenerative disc pathology, traumatic vertebral fractures, vertebral tumors, spinal stenosis, spondylolisthesis, pseudoarthrosis, and previous unsuccessful attempts at vertebral arthrodesis.

Any surgical decisions other than those recommended by the manufacturer are at the discretion and responsibility of the surgeon.

Do not use 4.5 mm diameter screws in the lumbar and lumbosacral spine, and do not couple 4.5 mm diameter screws with CoCr rods.



Features



TITANIUM ALLOY



ALLOY CO-CR
UPON REQUEST



ETO STERILE

#LESSISMORE

INSTRUMENTS



Clover has invested heavily in instrument design and care with the goal of creating ergonomic, functional, and compact instrumentation. Designed for the surgeon and his team.


TRAY 1



1 CEMENT NEEDLE ADAPTER	MSN-K0SS00200S	8 TISSUE DISSECTOR	MSN-T0SS00000S
2 DUAL LEAD CAP SCREWDRIVER	MSN-I1SS28525S	9 PIPE CLEANER	MSN-J2SS00000S
3 REAMING AWL	MSN-A1SS03507S	10 TULIP ALIGNMENT	MSN-K0SS00055S
4 SCREW TAB REMOVER	MSN-Z2SS00002S	11 RATCHETING T HANDLE	MSN-H1SS00100S
5 GUIDE WIRE TROCAR	MSN-K0NT00500S	12 RATCHETING HANDLE	MSN-H0SS00100S
6 TAP 4.5MM	MSN-J0S 500045S	13 DINAMOMETRIC T HANDLE 9NM	MSN-H1SS00090S
TAP 5.5MM	MSN-J0SS00055S	14 PIVOT 2	MSN-Z1SS00002S
TAP 6.5MM	MSN-J0SS00065S	15 PIVOT 1	MSN-Z1SS00001S
TAP 7.5MM	MSN-J0SS00075S	16 SUPERIOR RING	MSN-K2SS00015S
7 COUNTER TORQUE HANDLE	MSN-H2SS00000S	17 SLIDING RING	MSN-K2SS00025S

INSTRUMENTS

TRAY 2



18 MIS GAUGE	MSN-K3SS001385	24 DILATOR N.2	MSN-L2SS14409S
19 MIS ROD PUSHER	MSN-P1SS000015	25 DILATOR N.1	MSN-L0SS20002S
20 MIS TROCAR	MSN-A1SS02922S	26 DILATOR N.3	MSN-L2SS11514S
21 MIS COMPRESSOR	MSN-N0SS00002S	27 MIS CALIPER	MSN-M1SS00000S
22 ROD BENDER	MSN-F0SS00000S	28 DUAL LEAD MIS SCREWDRIVER	MSN-I2SS23520S
23 MIS ROD HOLDER	MSN-D1SS00000	29 CAP HOLDER	MSN-I0SS30025S

INSTRUMENTS

CEMENT NEEDLE ADAPTER

MSN-K0SS00200S



CAP SCREWDRIVER

MSN-I1SS28525S



DUAL LEAD MIS SCREWDRIVER

MSN-I2SS23520S



GUIDE WIRE TROCAR

MSN-KONT00500S



COUNTER TORQUE HANDLE

MSN-H2SS00000S



TAP 4.5MM

MSN-J0SS00045S



TAP 5.5MM

MSN-J0SS00055S



TAP 6.5MM

MSN-J0SS00065S



TAP 7.5MM

MSN-J0SS00075S



TAP 8.5MM

MSN-J0SS00085S



INSTRUMENTS

RATCHETING HANDLE

MSN-H0SS00100S



TULIP ALIGNMENT

MSN-K0SS00055S



RATCHETING T HANDLE

MSN-H1SS00100S



DINAMOMETRIC
T HANDLE 9Nm

MSN-H1SS00090S



MIS ROD HOLDER

MSN-D1SS00000



DILATOR N.1

MSN-L0SS20002S



DILATOR N.2

MSN-L2SS14409S



DILATOR N.3

MSN-L2SS11514S



MIS ROD PUSHER

MSN-P1SS00001S



MIS TROCAR

MSN-A1SS02922S



INSTRUMENTS

CAP HOLDER

MSN-I2SS24520S



REAMING AWL

MSN-A1SS13507S




MIS CALIPER

MSN-M1SS00000S




MIS COMPRESSOR

MSN-N0SS00002S



ROD BENDER

MSN-F0SS00000S



PIPE CLEANER

MSN-J2SS00000S



TISSUE DISSECTOR

MSN-T0SS00000S




MIS GAUGE

MSN-K3SS00138S



SLIDING RING

MSN-K2SS00025S



SUPERIOR RING

MSN-K2SS00015S





INSTRUMENTS

PIVOT 1

MSN-Z1SS00001S



PIVOT 2

MSN-Z1SS00002S




SCREW TAB REMOVER

MSN-Z2SS00002S



MIS SURGICAL TECHNIQUE

1 —



Preparation of the pedicle

After locating the access point, position the **trocar** and through radiographic control proceed to the inside of the pedicle.

Once the appropriate depth has been reached remove the core and handle, and insert the **guide wire** ensuring that it intercepts the vertebral body to ensure minimum anchorage to the vertebral body.

Then remove the trocar making sure not to move the guide wire.

Continue with the insertion of **dilator tube no. 1** and then **dilator tube no. 2** and **dilator tube no. 3**, which, thanks to the toothed termination allows anchorage to the articular process to prevent unwanted translation or displacement during manoeuvres.

MIS SURGICAL TECHNIQUE



Then remove **dilator tube No. 1**, checking that the guide wire remains in place, and should the operator wish to do so, prepare the screw housing with the **tapping tool**, which should be one size undersized in diameter from that of the screw you intend to place.

In case it is deemed appropriate, after removing **dilator tube No. 2** use the **reaming awl** to enable better positioning of the screw.

2 —



Screw assembly

Next proceed to the assembly of the screw.

Hook the **MIS screwdriver** onto the **ratcheting handle**.

Next assemble the screw to the **screwdriver** by inserting it inside the tulip and turn the ferrule clockwise to make the final tightening.

MIS SURGICAL TECHNIQUE

3 —



Screw insertion

Proceed to insert the screw inside the stalk following the **guide wire** and tighten.

Then slide out the **screwdriver**, turning counterclockwise the ferrule previously used to tighten the screw, remove the **guide wire** and **dilator tube No. 3**.

Insert the **sliding ring** along the tulip of the screw to prevent premature breakage of the screw, and if deemed appropriate, also insert the **superior ring** to prevent the screws from interfering with each other during subsequent operations.

Repeat the same operations for the insertion of subsequent screws.

Opt. —



Cementing


Connect the **cement needle adapter** by inserting its tip inside the screw, and then turn the guide ring clockwise to make the final tightening.

After checking that the two elements are firmly seated together, proceed to introduce the disposable bone filler to cement the screw.

Before injecting the cement, make sure that the bone filler has come to rest against the screw stem.

MIS SURGICAL TECHNIQUE

4 —



Rod insertion

After positioning the screws, perform the measurement of the distance between them using appropriate **MIS caliper** in order to choose the appropriate bar.

Please note that the gauge shows the actual measurement between the screw heads, so it is recommended to insert a bar with a length at least 5 mm longer than that shown on the **MIS caliper**.


If it is considered appropriate, use the **tissue dissector** to create the necessary space for inserting the bar.

Connect the bar to the **MIS rod holder** by inserting the hexagonal part of the bar into the appropriate hole on the gauge and secure it to the gauge by turning the locking crown clockwise.

If necessary, the bar can be bent and shaped using the **bar bending pliers**.

Then insert the bar by passing inside the tulips previously aligned according to a hypothetical sagittal plane.

5 —



Inserting the tightening nut

Using the **rod pusher** you can both check that the rod is seated correctly and push the rod down.

When you are certain that the bar has been properly seated in all screws, insert the nuts with the appropriate nut locators.

MIS SURGICAL TECHNIQUE




Once the nuts are in place, attach the **ratcheting handle** to the **cap screwdriver** and tighten the locknuts until the rod is pushed into the tulips.

Once the rod is secured to the screws, pull the rings out of the tulips.



6 —



Compression and distraction

Cover the tulips with **MIS gauge** and use **pivot 1** or **pivot 2** as the fulcrum of action in the center of the cannulas-after securing them to a handle.

To perform a compression, use the **MIS compressor** and compress below the **pivot**.

To perform a distraction, use the **MIS compressor** and compress above the **pivot**.



MIS SURGICAL TECHNIQUE

7 —




Final tightening

With the help of the **counter torque handle** connected to the **MIS gauge** make the final closure with the **9Nm dynamometric T-handle** assembled to the **cap screwdriver**.

Pull the **MIS gauge** out of the tulips and then unscrew the **rod holder** using the locking crown and then remove the gauge.

8 —



Tulip removal

Using the **screw tab removal** coupled with a **handle**, hook each individual tulip and then pry it loose from the screw cup.

To avoid displacements or small deformations of the system caused by the force required to break the tulips, a **counter torque** coupled with a handle can be used.

MONOAXIAL SCREW CANNULATED - FENESTRATED	D5.5 D6.5 D7.5 D8.5	FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L55
POLIAXIAL SCREW CANNULATED - FENESTRATED	D4.5 D5.5 D6.5 D7.5 D8.5 D9.5	FROM L25 TO L40 FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L90 FROM L30 TO L90 FROM L45 TO L90
TITANIUM PRECURVED ROD	D5.5	FROM L25 TO L100
TITANIUM ROD	D5.5	FROM L110 TO L500
CO-CR ROD	D5.5	FROM L110 TO L500
CLAMPING CAP FOR SCREWS AND HOOKS		
DOMINOES		
LONGITUDINAL CONNECTION ELEMENT		
MODULAR CROSSLINK		S - M - L - XL
ANGLED OFFSET		FROM L20 TO L60
OFFSET		FROM L20 TO L60
LUMBAR WIDE HOOK		FROM SIZE 7 TO SIZE 11
PEDICLE HOOK		FROM SIZE 5 TO SIZE 9
TORACHIC LAMINAR HOOK		FROM SIZE 5 TO SIZE 9
OLBIQUE HOOK - DX/SX		
OFFSET HOOK - DX/SX		
LUMBAR NARROW HOOK		FROM SIZE 7 TO SIZE 11

misano 

monza 

evo 

dixi 

recon-i 

rally 



Clover Orthopedics s.r.l.

Via Gadames n. 57/7, c.a.p. 20151 Milano

M. info@cloverorthopedics.com

W. cloverorthopedics.com

T. +39 02 457 902 31

F. +39 02 457 902 66

CE
0426

 **ITALCERT**
SISTEMA DI GESTIONE QUALITÀ
CERTIFICATO
UNI CEI EN ISO 13485:2016
UNI EN ISO 9001:2015