

misano

CLOVER





Complete set of sizes for open, mis and deformity procedures

Self-tapping screw for insertion even without tapping

Great mechanical properties

Double thread for improved sealing

Sophisticated yet compact instrumentation

Cannulated screw for injection of cement or bone substitute


Sterile packaging

Stabilization system thoracolumbar vertebral

Misano is the first stabilization system thoracolumbosacral, totally made of titanium, which allows the use of a single pedicle screw for the treatment of all degenerative, traumatic, and deformity pathologies. The self-tapping screw can be inserted without the aid of the tap, has a specific two-principle thread at the distal level that ensures rapid implantation, which then changes to a four-principle thread also ensuring a high pedicle seal. The screws are also always cannulated to allow insertion with the aid of a guide wire and fenestrated to allow cementing.

The Misano thoracolumbosacral stabilization system from Clover Orthopedics consists of single- and polyaxial pedicle screws, clamping nuts, straight and precurved titanium and cobalt-chrome bars, a wide range of hooks, offsets, connectors, and cross-links, and is used to facilitate a solid arthrodesis of the treated portion of the spine.





Appropriately used, the Misano thoracolumbar-sacral stabilization system from Clover Orthopedics is indicated to promote the development of solid thoracic, lumbar and sacral arthrodesis.

It is recommended in cases of spinal deformity, degenerative disc disease, traumatic vertebral fractures, vertebral tumors, spinal stenosis, spondylolisthesis, pseudoarthrosis, and previous unsuccessful attempts at vertebral arthrodesis.

Any surgical decisions other than those recommended by the manufacturer are at the discretion and responsibility of the surgeon.

Do not use 4.5 mm diameter screws in the lumbar and lumbosacral spine, and do not couple 4.5 mm diameter screws with CoCr rods.

Features



TITANIUM ALLOY



ALLOY CO-CR
UPON REQUEST



ETO STERILE




#LESSISMORE

INSTRUMENTS



Clover has invested heavily in instrument design and care with the goal of creating ergonomic, functional, and compact instrumentation. Designed for the surgeon and his team.


CASE 1 | TRAY 1



1	THORACIC AWL	MSN-A0SS02510S	8	T-HANDLE	MSN-H1SS00000S
2	LUMBAR AWL	MSN-A0SS03208S	9	RATCHETING T-HANDLE	MSN-H1SS00100S
3	NARROW GEAR SHIFT PROBE CURVED	MSN-B1SS00000S	10	ROD TEMPLATE L200 D5.5	MSN-M0AL00000S
4	NARROW GEAR SHIFT PROBE STRAIGHT	MSN-B0SS00000S	11	SOUNDING 2MM	MSN-C0SS00000S
5	GEAR SHIFT PROBE STRAIGHT	MSN-B0SS00001S	12	DOUBLE SOUNDING 2MM	MSN-C1SS00000S
6	GEAR SHIFT PROBE CURVED	MSN-B1SS00001S	13	KIRCHNER WIRES	MSN-KONT00500S
7	RATCHETING HANDLE	MSN-H0SS00100S			

INSTRUMENTS


CASE 1 | TRAY 2



14	POLYAXIAL SCREWDRIVER	MSN-I1SS25120S	22	FREEHAND POLYAXIAL SCREWDRIVER	MSN-I1SS28520S
15	POLYAXIAL CANNULATED SCREWDRIVER	MSN-I2SS23920S	23	TAP 4.5MM	MSN-J0SS00045S
16	CASPAR	MSN-N0SS00001S		TAP 5.5MM	MSN-J0SS00055S
17	CAP HOLDER	MSN-I0SS30025S		TAP 6.5MM	MSN-J0SS00065S
18	DINAMOMETRIC T-HANDLE 9Nm	MSN-H1SS00090S		TAP 7.5MM	MSN-J0SS00075S
19	ROD HOLDER	MSN-D0SS00001S	24	CEMENT NEEDLE ADAPTER	MSN-K0SS00200S
20	CAP SCREWDRIVER	MSN-I1SS28525S	25	COUNTER TORQUE HANDLE	MSN-Q0SS00000S
21	TULIP ALIGNMENT	MSN-K0SS00227S	26	LOCKSCREW GUIDE TOWER	MSN-Q0SS00002S

INSTRUMENTS

CASE 2 | TRAY 1



27	CAP HOLDER LONG	MSN-I05S380255	32	ROD BENDER	MSN-F05S000005
28	HAEXAGONAL WRENCH	MSN-Q05S000015	33	ROD PUSHER	MSN-P05S000005
29	PERSUADER	MSN-P05S000055	34	DISTRACTOR	MSN-N15S000005
30	ROCKER	MSN-P05S000035	35	COMPRESSOR	MSN-N05S000005
31	ROD GRIPPER	MSN-E15S000005			

INSTRUMENTS

THORACIC AWL

MSN-A0SS02510S



LUMBAR AWL

MSN-A0SS03208S



NARROW GEAR SHIFT
PROBE CURVED

MSN-B1SS00000S



NARROW GEAR SHIFT
PROBED STRAIGHT

MSN-B0SS00000S



GEAR SHIFT
PROBE CURVED

MSN-B1SS00001S



GEAR SHIFT
PROBE STRAIGHT

MSN-B0SS00001S



RATCHETING HANDLE

MSN-H0SS00100S



RATCHETING T HANDLE

MSN-H1SS00100S



T HANDLE

MSN-H1SS00000S



ROD TEMPLATE L200 D5.5

MSN-M0AL00000S



INSTRUMENTS

SOUNDING 2MM

MSN-C0SS000005



SOUNDING 2MM DOUBLE

MSN-C1SS000005



POLYAXIAL SCREWDRIVER

MSN-I1SS251205



POLYAXIAL CANNULATED SCREWDRIVER

MSN-I2SS239205



FREEHAND POLYAXIAL SCREWDRIVER

MSN-I1SS285205



CAP HOLDER

MSN-I0SS300255




CAP SCREWDRIVER

MSN-I1SS285255



CASPAR

MSN-N0SS000015



ROD HOLDER

MSN-D0SS000015



DINAMOMETRIC
T HANDLE 9Nm

MSN-H1SS000905



INSTRUMENTS

TULIP ALIGNMENT MSN-K0SS00227S



TAP 4.5MM MSN-J0SS00045S



TAP 5.5MM MSN-J0SS00055S



TAP 6.5MM MSN-J0SS00065S



TAP 7.5MM MSN-J0SS00075S



KIRCHNER WIRES MSN-KONT00500S



COUNTER-TORQUE HANDLE MSN-Q2SS00000S



LOCKSCREW GUIDE
DEROTATION TOWER MSN-Q0SS00002S



ROD PUSHER MSN-P0SS00000S




HEAXGONAL WRENCH MSN-Q0SS00001S



INSTRUMENTS


ROD BENDER

MSN-F0SS000005




ROD GRIPPER

MSN-E1SS000005




ROCKER FORCEPS

MSN-P0SS000035



PERSUADER

MSN-P0SS000055




CEMENT NEEDLE ADAPTER

MSN-K0SS002005




PARALLEL DISTRACTOR

MSN-N1SS000005




PARALLEL COMPRESSOR

MSN-N0SS000005



OPEN SURGICAL TECHNIQUE

1 —



Preparation of the pedicle screw

After determining the entry point proceed via **thoracic or lumbar awl** through the pedicle.

Subsequently, open the pedicle canal with the **curved or straight gear shift**.


Use the **straight or curved sounding** to assess the integrity of the pedicle walls.

Although Misano pedicle screws are self-tapping, **taps** can be used to prepare the screw housing. Tapping undersized by one size relative to the selected screw size is recommended. The taps - available in diameters 4.5 mm, 5.5 mm, 6.5 mm, 7.5 mm and 8.5 mm - must be assembled to the **handle**. The thread length is 30 mm.

All of the above tools are laser-marked at 10 mm intervals to display depth and aid in selecting the appropriate screw length.

OPEN SURGICAL TECHNIQUE

2 —




Inserting the screw

Connect the **polyaxial screwdriver** to the **ratcheting handle**:

- connect the **screwdriver** to the screw by inserting its tip inside the tulip until the thread reaches the head of the screw and then turn the screwdriver nut clockwise to make the final tightening;
- proceed to implant the screw.

Once the desired insertion depth is reached, the screwdriver is disengaged from the screw by turning the screwdriver nut counterclockwise.

3 —



Adjustment of the screw head

If necessary, you can use the **freehand screwdriver** by placing the tip of the screwdriver directly in contact with the stem of the screw. Tighten or unscrew the screw to the desired depth.

The head of the polyaxial screw can be aligned using the **tulip alignment**.

OPEN SURGICAL TECHNIQUE

Opt. —



Cementing

Connect the **cement needle adapter** by inserting its tip inside the screw, and then turn the guide ring clockwise to make the final tightening.

After checking that the two elements are firmly seated together, proceed to introduce the disposable bone filler to cement the screw.

Before injecting the cement, make sure that the bone filler has come to rest against the screw stem.

4 —



Rod preparation

Once all screws are in place use the **rod template** to determine the length of the final bar.

If necessary to bend and shape the bars you can use the **rod bender**.



5 —




Positioning the rod

Place the previously selected bar inside the screw head with the **rod holder**.

OPEN SURGICAL TECHNIQUE

6A —



Positioning of the clamping nut — A

To ensure proper positioning of the clamping nut on the screw, the **lock screw guide derotation tower** can be used. The latter has a slightly larger distal end that is positioned over the head of the screw. At the same time, by exerting downward pressure, you properly seat the bar inside the screw. You can then proceed to position the clamping nut with the **cap holder**.

If necessary, the **rod pusher** can be used.



6B —




Positioning of the clamping nut — B

As an alternative to the **rod pusher**, the **rocker** can be used. Grasp the waist head with the rocker, subsequently deflect the rocker downward until the bar is fully inserted into the screw head. Then insert a clamping nut using the **cap holder**.



OPEN SURGICAL TECHNIQUE

6C —



Positioning of the clamping nut — C


Alternatively, it is possible to use the **axial persuader with folding handle**, which allows, with greater force, the bar to be brought closer to the head of the screw and the tightening nut then positioned.

With the handles fully open, engage the persuader to the tulip of the screw. By then tightening the handles, the mechanism lowers the bar until it is inserted into the head of the screw.

Next, introduce the clamping nut inside the persuader using the long clamping **cap holder**.

Then release the persuader from the screw tulip by releasing the handles.

7 —




Rod rotation

If it is necessary to rotate the bar before final clamping, use the **rod gripper**.

Alternatively, the hexagonal end of the bar allows the additional use of the **hexagonal wrench**.

OPEN SURGICAL TECHNIQUE

8 —



Compression and Distraction

Once the bar has been fully inserted into all polyaxial screw heads and the clamping nuts are in place, distraction or compression can be performed using the **distractor** or **compressor**.

9 —



Final tightening

Proceed to the final tightening of the system.

Final tightening is carried out using the **countertorque** system, consisting of the **countertorque handle** connected to the **lock screw guide derotation tower**, and the **cap screwdriver** previously assembled to the **9Nm dynamometric T-handle**.

MONOAXIAL SCREW CANNULATED - FENESTRATED	D5.5 D6.5 D7.5 D8.5	FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L55
POLIAXIAL SCREW CANNULATED - FENESTRATED	D4.5 D5.5 D6.5 D7.5 D8.5 D9.5	FROM L25 TO L40 FROM L30 TO L55 FROM L30 TO L55 FROM L30 TO L90 FROM L30 TO L90 FROM L45 TO L90
TITANIUM PRECURVED ROD	D5.5	FROM L25 TO L100
TITANIUM ROD	D5.5	FROM L110 TO L500
CO-CR ROD	D5.5	FROM L110 TO L500
CLAMPING CAP FOR SCREWS AND HOOKS		
DOMINOES		
LONGITUDINAL CONNECTION ELEMENT		
MODULAR CROSSLINK		S - M - L - XL
ANGLED OFFSET		FROM L20 TO L60
OFFSET		FROM L20 TO L60
LUMBAR WIDE HOOK		FROM SIZE 7 TO SIZE 11
PEDICLE HOOK		FROM SIZE 5 TO SIZE 9
TORACHIC LAMINAR HOOK		FROM SIZE 5 TO SIZE 9
OLBIQUE HOOK - DX/SX		
OFFSET HOOK - DX/SX		
LUMBAR NARROW HOOK		FROM SIZE 7 TO SIZE 11

misano 

monza 

evo 

dixi 

recon-i 

rally 



Clover Orthopedics s.r.l.

Via Gadames n. 57/7, c.a.p. 20151 Milano

M. info@cloverorthopedics.com

W. cloverorthopedics.com

T. +39 02 457 902 31

F. +39 02 457 902 66

CE
0426

 SISTEMA DI GESTIONE QUALITÀ
CERTIFICATO
ITALCERT
UNI CEI EN ISO 13485:2016
UNI EN ISO 9001:2015