


● CLOVER

monza
CERVICAL





Possibility of stand-alone or stabilized approach

Specific tool for easy and safe removal of the anchors

Zero-profile cage with stabilizing blades that preserve bone tissue compared to the use of screws

Trabecular structure to facilitate the process of fusion between the vertebral bodies

Cervical cage in trabecular titanium

Monza Cervical is an innovative cervical cage that, thanks to its titanium trabecular structure made with the latest 3D printing techniques, provides immediate and secure mechanical stability and certain osseointegration to all types of implants.

The cervical cage was developed to cover all possible sizes and degrees of lordosis and allows placement with or without the use of "anti-migration" ratchets.

The instrumentation is basic but remarkably effective and integrated with the cage when inserted and positioned.





Features



TRABECULAR TITANIUM



3D PRINTING TECH



ETO STERILE


Instrumentary

Clover has invested heavily in instrument design and care with the goal of creating ergonomic, functional, and compact instrumentation. Designed for the surgeon and his team.



#LESSISMORE

TRAY 1



1	SCRAPER	MNC-A0SS000005	6	REMOVAL HOOK	MNC-L0SS000005
2	CURETTE	MNC-B0SS000005	7	IMPACTOR	MNC-G0SS000005
3	TEARDROP CURETTE	MNC-B0SS000015	8	HEX DRIVER	MNC-H0SS000005
4	TRIAL 12X14 H5-H6	MNC-C0SS005065	9	LOCKING SYSTEM HOLDER	MNC-H1SS000005
	TRIAL 12X14 H7-H8	MNC-C0SS007085	10	POSTERIOR CAP	MNC-I0SS000005
	TRIAL 12X14 H9-H10	MNC-C0SS009105	11	RATCHET HOLDER	MNC-D3SS000005
	TRIAL 14X16 H5-H6	MNC-C2SS005065	12	HOLDER STAND ALONE	MCN-D2SS000005
	TRIAL 14X16 H7-H8	MNC-C2SS007085	13	IMPLANT HOLDER WITHOUT STOP	MNC-D1SS000005
	TRIAL 14X16 H9-H10	MNC-C2SS009105	14	HOLDER WITH SAFE STOP	MNC-D0SS000005
5	SLIDE HAMMER	MNC-E0SS000005			

INSTRUMENTS

SCRAPER

MNC-A0SS000005



CURETTE

MNC-B0SS000005



TEARDROP CURETTE

MNC-B0SS000015



TRIAL 12X14 H5-H6

MNC-C0SS005065



TRIAL 12X14 H7-H8

MNC-C0SS007085



TRIAL 12X14 H9-H10

MNC-C0SS009105



TRIAL 14X16 H5-H6

MNC-C2SS005065



TRIAL 14X16 H7-H8

MNC-C2SS007085



TRIAL 14X16 H9-H10

MNC-C2SS009105



REMOVAL HOOK

MNC-L0SS000005



INSTRUMENTS

IMPACTOR MNC-G0SS000005



HEX DRIVER MNC-H0SS000005



LOCKING SYSTEM HOLDER MNC-H1SS000005



RATCHET HOLDER MNC-D3SS000005



HOLDER STAND ALONE MCN-D2SS000005



IMPLANT HOLDER WITHOUT STOP MNC-D1SS000005



IMPLANT HOLDER WITH SAFE STOP MNC-D0SS000005



POSTERIOR CAP MNC-I0SS000005




SLIDE HAMMER MNC-E0SS000005



SURGICAL TECHNIQUE

1 —




Disc space preparation

If necessary, use the distractor to achieve adequate access to the disc space. The pins should be placed approximately 7 mm from both vertebral plates to avoid contact between the caspar pins and the Monza anchor flap during insertion of the anchor flap.

If the patient's anatomy does not allow a space of 7 mm, remove the caspar pins before inserting the anchor fin to eliminate any risk of obstruction.

Remove the disc and perform any removal of bone and tissue with special instruments found within the instrumentarium: **curettes** and **scrapers**.

2 —



Cage selection

Select the appropriate cage using the **trial** provided.

NOTE: Evidence and cages have an anatomical profile.

3 —






Preparation and insertion of the standalone cage

Screw the cage onto the **implant holder** and insert it into the intervertebral disc space. Then, using fluoroscopy, check the correct positioning of the cage, and if the position is correct, remove the cage holder.

SURGICAL TECHNIQUE

4 —



Preparation and insertion of the cage with ratchet

Screw the cage onto the **implant holder** by turning the knob on the end of the holder. Then with the help of the protective ring, insert the cage into the intervertebral space and, using fluoroscopy, check the correct positioning of the cage.

Insert the two ratchet tabs into the groove of the specific implant holder; the direction of the tabs is defined on the top of the cage holder by laser marking arrows.

Using the **impactor**, insert the ratchets inside the vertebral bodies. It is advisable to insert one ratchet at a time and verify with fluoroscopy the correct positioning of the fin. When the impactor makes contact with the top of the **implant holder**, it means that the ratchet is positioned correctly and totally inside the vertebral body.

Remove the **implant holder** by turning the knob on the end of the holder.

With the nut locator, screw the nut into the hole provided to prevent the lugs from sticking out of the cage. With the nut positioner screwdriver, a greater closing force can be exerted.

CERVICAL LORDOTIC CAGE - 6°	12X14 14X16	FROM H5 TO H10 FROM H5 TO H10
RATCHET LOCKING SYSTEM		H5-H6 H7-H8 H9-H10

misano 

monza 

evo 

dixi 

recon-i 

rally 



Clover Orthopedics s.r.l.
Via Gadames n. 57/7, c.a.p. 20151 Milano

M. info@cloverorthopedics.com
W. cloverorthopedics.com

T. +39 02 457 902 31
F. +39 02 457 902 66

CE
0426

 **ITALCERT**
SISTEMA DI GESTIONE QUALITÀ
CERTIFICATO
UNI CEI EN ISO 13485:2016
UNI EN ISO 9001:2015