


● CLOVER



monza  
LUMBAR



ZERO PROFILE CAGE, USE OF STABILIZING BLADES, WHICH PRESERVE BONE TISSUE COMPARED WITH THE USE OF SCREWS

SLOT FOR THE INSERTION OF AN ADDITIONAL STABILIZATION

TRABECULAR STRUCTURE TO FACILITATE THE FUSION PROCESS BETWEEN THE VERTEBRAL BODIES

DEDICATED INSTRUMENT FOR INSERTION CAGE WITH A GUIDE WIRE

## Complete system of cages inter-somatic thoracolumbar

Monza is a complete system of interbody cages that, thanks to the titanium trabecular structure, made with the latest 3D printing techniques, guarantees immediate mechanical stability and excellent osseointegration.

The Plif, Tlif and Llif lumbar cages come in different sizes and degrees of lordosis ensuring a wide modularity of the implant.

The design of the implants leaves ample room for the insertion of bone tissue .

The instrumentation is ergonomic, essential and effective.






TLIF



TLIF RACE



PLIF

## Features



TITANIO TRABECOLARE



3D PRINTING TECH



ETO STERILE


## Instrumentary

Clover has invested heavily in instrument design and care with the goal of creating ergonomic, functional, and compact instrumentation. Designed for the surgeon and his team.



## INSTRUMENTS


### TRAY 1



<b>1</b>	HOLDER	MNZ-B0SS00000S	<b>7</b>	ANGLED TEARDROP CURETTE	MNZ-C1SS00001S
<b>2</b>	CHISEL	MNZ-P0SS00004S	<b>8</b>	ROTATING TLIF HOLDER	MNZ-B0SS00001S
<b>3</b>	CURETTE	MNZ-C0SS00000S	<b>9</b>	CANNULATED T-HANDLE	MNZ-N1SS00000S
<b>4</b>	TEARDROP CURETTE	MNZ-C0SS00001S	<b>10</b>	MOULD FOR BONE GRAFT	MNZ-F2SS00000S
<b>5</b>	STRAIGHT SCRAPER	MNZ-D0SS00000S	<b>11</b>	IMPACTOR FOR BONE GRAFT	MNZ-G0SS00000S
<b>6</b>	ANGLED SCRAPER	MNZ-D1SS00000S			

## INSTRUMENTS

### TRAY 2



12	EXTRACTOR	MNZ-L0SS000005	16	SLIDE HAMMER	MNZ-I0SS000005
13	STARTER 7MM	MNZ-A0SS000005	17	LATERAL IMPACTOR	MNZ-M0SS000015
14	SPREADER / TRIAL PL-STL 8MM	MNZ-H0SS000085	18	POSTERIOR IMPACTOR	MNZ-M0SS000005
	SPREADER / TRIAL PL-STL 10MM	MNZ-H0SS000105	19	PL-TL INSERT	MNZ-P0SS000055
	SPREADER / TRIAL PL-STL 12MM	MNZ-H0SS000125	20	TRIAL TL 7MM	MNZ-E0SS000075
	SPREADER / TRIAL PL-STL 14MM	MNZ-H0SS000145		TRIAL TL 8MM	MNZ-E0SS000085
15	SHAVER 7MM	MNZ-O0SS000075		TRIAL TL 9MM	MNZ-E0SS000095
	SHAVER 8MM	MNZ-O0SS000085		TRIAL TL 10MM	MNZ-E0SS000105
	SHAVER 9MM	MNZ-O0SS000095		TRIAL TL 11MM	MNZ-E0SS000115
	SHAVER 10MM	MNZ-O0SS000105		TRIAL TL 12MM	MNZ-E0SS000125
	SHAVER 11MM	MNZ-O0SS000115		TRIAL TL 13MM	MNZ-E0SS000135
	SHAVER 12MM	MNZ-O0SS000125			
	SHAVER 13MM	MNZ-O0SS000135			
	SHAVER 14MM	MNZ-O0SS000145			



## INSTRUMENTS

HANDLE MNZ-N05S000005



T-HANDLE MNZ-N15S000005



CHISEL MNZ-P05S000045



ANGLED TEARDROP CURETTE MNZ-C15S000015



TEARDROP CURETTE MNZ-C05S000015



CURETTE MNZ-COSS000005



STRAIGHT SCRAPER MNZ-D05S000005



ANGLED SCRAPER MNZ-D15S000005



HOLDER MNZ-B05S000005



ROTATING TLIF HOLDER MNZ-B05S000015



## INSTRUMENTS

IMPACTOR FOR BONE GRAFT MNZ-G0SS000005



MOULD FOR BONE GRAFT  
TLIF CAGES MNZ-F1SS000005



POSTERIOR IMPACTOR MNZ-M0SS000005



LATERAL IMPACTOR MNZ-M0SS000015



EXTRACTOR PLIF-STLIF MNZ-L0SS000005



SHAVER 7MM MNZ-O0SS000075



SHAVER 8MM MNZ-O0SS000085



SHAVER 9MM MNZ-O0SS000095



SHAVER 10MM MNZ-O0SS000105



SHAVER 11MM MNZ-O0SS000115



## INSTRUMENTS

SHAVER 12MM

MNZ-00SS00012S



SHAVER 13MM

MNZ-00SS00013S



SHAVER 14MM

MNZ-00SS00014S



TRIAL TL 7MM

MNZ-E0SS00007S



TRIAL TL 8MM

MNZ-E0SS00008S



TRIAL TL 9MM

MNZ-E0SS00009S



TRIAL TL 10MM

MNZ-E0SS00010S



TRIAL TL 11MM

MNZ-E0SS00011S



TRIAL TL 12MM

MNZ-E0SS00012S



TRIAL TL 13MM

MNZ-E0SS00013S





## INSTRUMENTS

TRIAL PL/S-TL 7MM

MNZ-A0SS000005



TRIAL PL/S-TL 8MM

MNZ-H0SS000085



TRIAL PL/S-TL 10MM

MNZ-H0SS000105



TRIAL PL/S-TL 12MM

MNZ-H0SS000125



TRIAL PL/S-TL 14MM

MNZ-H0SS000145



PLIF-TLIF INSERT

MNZ-P0SS000055






SLIDE HAMMER

MNZ-I0SS000005



## PLIF / S-TLIF SURGICAL TECHNIQUE

1 —




### Disc space preparation

After identifying the implant site perform resection of the affected anatomic portion via posterior access.

Once the disc is reached, proceed with discectomy using sized **shavers**, **curettes** and **scrapers** in addition to the appropriate general instruments. The size of the last shaver used may provide a useful indication for the subsequent choice of cage trial.

2 —



### Selection of the cage

Then proceed to insert the **PL/S-TL trial** inside the disc, starting with the smallest size until the desired height is reached.

NOTE: There are grooves visible on the end of the trial under scopy to determine the length of the cage to be selected.


3 —



### Preparing the cages

Mount the cages on the **holder**, ensuring assembly by turning the ring nut on the handle clockwise.

## PLIF / S-TLIF SURGICAL TECHNIQUE



Place the cages in the **mould for bone graft PLIF/S-TLIF cages** and introduce the bone graft inside the implant hole using the appropriate **impactor for bone graft**.

4 —




## Insertion of cages


Then insert the cages inside the disc space to the desired depth. If you are not satisfied with the position achieved, you can connect the impactor directly to the positioner and remove the cage.

Once the desired depth and position has been reached, unscrew the ferrule and remove the positioner.

With the positioner removed, the PL/S-TL definitive beater can be used to slightly advance the cage within the disc space.



5 —




## Cage removal

Screw the **extractor** inside the cage hole. Connect the instrument to the **slide hammer** and proceed to extract the cage.

## TLIF SURGICAL TECHNIQUE

1 —




### Disc space preparation

After identifying the implant site perform resection of the affected anatomic portion via transforaminal access.

Once the disc is reached proceed with discectomy using, in addition to the appropriate general instruments, sized **shavers**, **curettes**, and **scrapers**. The size of the last shaver used may provide a useful indication for the subsequent choice of cage trial.

2 —




### Selection of the cage

Then proceed to insert the **TL trial** inside the disc, starting with the smallest size until the desired height is reached.

NOTE: The length of the TL test measures 28 mm.

3 —



### Preparing the cages

Mount the cages on the **holder** and turn the ring nut on the handle clockwise.


Next, in order to secure the cage to the **holder**, slide the center body of the **holder** into the cage and rotate counterclockwise.

## TLIF SURGICAL TECHNIQUE



Place the cages in the **mould for bone graft PLIF/S-TLIF cages** and introduce the bone graft inside the implant hole using the appropriate **impactor for bone graft**.

4 —



## Insertion of cages

Then insert the cages inside the disc space to the desired depth. If you are not satisfied with the position achieved, you can connect the impactor directly to the positioner and remove the cage.

Once the desired depth and position has been reached, unscrew the ferrule and remove the positioner.

With the positioner removed, the TL definitive beater can be used to gradually rotate and advance the cage under radiographic control to its final position.

PLIF CAGE	10X22	FROM H7 TO H14
PLIF CAGE LORDOTIC	10X26	FROM H7 TO H14
TLIF CAGE	12X28 12X32	FROM H7 TO H14 FROM H7 TO H14
S-TLIF CAGE	10X30	FROM H7 TO H14
LLIF CAGE	18X40 18X45 18X50 18X55	FROM H7 TO H12 FROM H7 TO H12 FROM H7 TO H12 FROM H7 TO H12
LLIF CAGE LORDOTIC- 10°	18X40 18X45 18X50 18X55	FROM H7 TO H14 FROM H7 TO H14 FROM H7 TO H14 FROM H7 TO H14

misano 

monza 

evo 

dixi 

recon-i 

rally 



Via Gadames n. 57/7, c.a.p. 20151 Milano

E. [info@cloverorthopedics.com](mailto:info@cloverorthopedics.com)  
W. [cloverorthopedics.com](http://cloverorthopedics.com)

T. +39 02 457 902 31  
F. +39 02 457 902 66

CE  
0426

 SISTEMA DI GESTIONE QUALITÀ  
CERTIFICATO  
**ITALCERT**  
UNI CEI EN ISO 13485:2016  
UNI EN ISO 9001:2015